

BAUERFEIND.

BY PEOPLE FOR PEOPLE. FOR A BETTER QUALITY OF LIFE.



Medical Device Establishment Licence Liscence No 4678 7017006 CANADA INC. O/A GALIEN

Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products and treatment concepts, we enable people to be more active and thus make an important contribution toward maintaining and restoring health.

MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE

Bauerfeind combines decades of expertise in technology and manufacturing with innovation: the tradition of a medium-sized, family-run business with a global mindset. Bauerfeind's specialists have worked very closely with physicians, scientists, and developers right from the start. The result: patented, award-winning products and services.

From the thread to the completed product, our compression stockings are developed and manufactured in Germany.

EFFECTIVE VARIETY

Our VenoTrain compression stockings cover all uses/indications in modern compression therapy, offering medically proven effectiveness and unique material properties.

Thanks to the wide range of colors available, the compression products blend seamlessly into their users' everyday lives and can be combined with various fashion styles. For prevention or the treatment of lymphedema: the indication-specific solutions include a wide variety of circular knit and flat knit designs that are available as custom-made or standard products.

QUALITY WITHOUT COMPROMISE

Bauerfeind products don't simply meet standards – they fulfill the most stringent requirements in terms of quality, effectiveness, and wearing comfort. Regular inspection procedures ensure that our products earn the seals and marks of successful testing.

Canadian distributor Galien

P Office toll free 1-877-629-9889 **F** Fax toll free 1-877-629-9895

E info@bauerfeind.ca





















CONTENTS

OVERVIEWS

- 06 The three pillars of compression therapy
- 07 A systematic approach to compression

PILLAR 01



09 PHLEBOLOGICAL
TREATMENT
STANDARD CASES

- 11 Compression Sock TRAINING
- 13 Compression Sock
 PERFORMANCE
- 15 **VenoTrain® discrétion**
- 17 VenoTrain® look
- 19 VenoTrain® micro
- 21 Design and Batik
- 23 VenoTrain® business

PILLAR 02



25 PHLEBOLOGICAL TREATMENT SPECIAL CASES

- 27 VenoTrain® soft
- 29 VenoTrain® soft arm
- 31 VenoTrain® soft S
- 33 VenoTrain® impulse
- 35 VenoTrain® ulcertec

PILLAR 03



37 LYMPHOLOGICAL TREATMENT

ACCESSORIES

- 39 Care accessories
- 40 VenoTrain®glider
- 41 VenoTrain® glider plus

THE THREE PILLARS OF COMPRESSION THERAPY AN ALL-IN-ONE SOLUTION

Bauerfeind products meet particularly high standards, allowing you to provide your customers with the best possible treatment. VenoTrain compression stockings offer indication-specific solutions for all areas of treatment in modern compression therapy: whether for prevention, for phlebological standard or special cases, or for conservative treatment of lymphedema or lipedema. The circular knit or flat knit compression stockings offer medically proven effectiveness, unique material properties, and a patient-oriented modern design.



TREATMENT IN PHLEBOLOGICAL STANDARD CASES

The phlebological standard case refers to **mild to moderate venous disorders**. It also includes situations in which there is a particular risk for the development of vein problems, such as during pregnancy or traveling. VenoTrain compression stockings improve venous activity and combat venous leg problems.



02 TREATMENT IN PHLEBOLOGICAL SPECIAL CASES

Some patients with vein problems have **special requirements** that call for compression stockings with unique modes of action and properties. These requirements can be very diverse. That is why the VenoTrain range offers compression stockings especially for patients with restricted mobility, dry skin, particularly weak veins, or open wounds




03 LYMPHOLOGICAL TREATMENT

Lymphedema or lipedema patients are generally dependent on compression therapy for their entire lives. The effectiveness and wearing comfort of the stockings is therefore particularly important, significantly impacting the patient's quality of life.

VenoTrain flat knit and circular knit products provide the required therapeutic compression, while offering a unique wearing comfort, and help the patient live an active life.

Detailed information about products for lymphedema and lipedema can be found in the separate Lymphology Product Catalog.

SYSTEMATIC APPROACH TO COMPRESSION



^{*} Partsch, H.: Klassifikation von chronischen Venenerkrankungen heute, Vasomed 1/2003

^{**}Rabe, E. et al.: Bonner Venenstudie der DGP, Phlebologie 1 / 2003



TREATMENT IN PHLEBOLOGICAL STANDARD CASES FOR MILD TO MODERATE VEIN PROBLEMS

The phlebological standard case refers to **mild to moderate venous disorders**. It also includes situations in which there is a particular risk for developing vein problems, such as during pregnancy or traveling. VenoTrain compression stockings have been designed precisely to improve venous activity and combat venous leg problems.



Color	Product Code	Length	Size
Coal/Rivera	291800530000	2 = Short	1 = Small
Silver/Rivera	291800531000	3 = Long	2 = Medium
Coal/Polar	291800532000		3 = Large
Silver/Polar	291800533000		4 = Xlarge
Coal/Pink	291800197000		
Silver/Pink	291800196000		

COMPRESSION SOCK **TRAINING**

MUSCLE TUNING MADE EASY!

Bauerfeind's Compression Sock Training supports your calf muscle and lower leg, and stimulates the parts of the muscles that are under particular stress. In addition to providing special benefits, the socks help to prevent premature fatigue and muscle cramping. The result: a lower risk of injury and more energy throughout your activities.

DETAILS

TOE closed **MATERIAL CONTENT** 93 % polyamide, 7 % elastane **CARE INSTRUCTIONS**



FEATURES

- Set pressure profile consistent therapeutic gradient (Ccl 1 = 20-30 mmHg) for improved circulation
- Lycra® sports material gentle on the skin and breathable
- Taping zone more stability for the ankle thanks to a patented design
- Precision fit anatomically correct for left and right leg
- Muscle tuning zone for reduced muscle vibrations and thus less fatigue
- Protective zone for less pressure and rubbing on the Achilles tendon
- Shock-absorbing cushioning for greater comfort and protection

AVAILABLE IN



	CIRCUMFE	RENCE	S	M	L	XL
Length ID						
short up to 41 cm (16")	cC	cm	31 – 40	34-44	37-44	40-51
long over 41 cm (16")	Calf	inch	12-15	13-17	14.5-18.5	15.5-20
Type Knee high						
	сВ	cm	20-22	22-25	25-28	28-31
	Ankle	inch	8 – 9	9 – 10	10-11	11-12

C0 C1 C2 C3 C4 C5 C6



	Ccl	Color	Тое	Circumference	Length	Size
292 800 = Knee high	1 = 20-30 mmHg	10 = White	00 = Closed toe	0 = Normal	2 = Short	1 = Small
		17 = Black			3 = Long	2 = Medium
		34 = Rivera				3 = Large
		97 = Pink				4 = Xlarge
		46 = Wild Berry				

COMPRESSION SOCK **PERFORMANCE**

RUNNING MADE EASY!

The Bauerfeind Compression Sock Performance exerts noticeable compression on the leg that complies with medical standards and effectively improves circulation. This ensures an optimal blood supply to the muscles and helps them perform for longer, particularly in endurance disciplines such as cycling or marathons.

The result: More energy reserves for long distances and quicker regeneration.

DETAILS

TOE closed
MATERIAL CONTENT
75% polyamide, 25% elastane
CARE INSTRUCTIONS



FEATURES

- Set pressure profile consistent therapeutic gradient (Ccl 1 = 20-30 mmHg) for improved circulation
- Calf stabilization reduces typical muscle vibrations and prevents premature fatigue
- Breathable microfiber provides optimum thermoregulation and is hard-wearing and comfortably light to wear
- **Elegant design** understated for professionals and goes perfectly with any style of clothing
- Non-skid foot design comfortable protection from pressure points and during increased strain

AVAILABLE IN



	CIRCUMFER	ENCE	S	M	L	XL
Length lD						
hort up to 41 cm (16")	cC	cm	31 – 40	34-44	37 – 44	40-51
ong over 41 cm (16")	Calf	inch	12-15	13-17	14.5 – 18.5	15.5-20
ype Knee high						
	сВ	cm	20 – 22	22-25	25 – 28	28-31
	Ankle	inch	8-9	9-10	10-11	11-12

C0 C1 C2 C3 C4 C5



	Ccl	Color	Toe	Circumference	Length	Size
216 8002 = Knee high	1 = 20-30 mmHg	02 = Clair	02 = Closed toe	2 = Normal	2 = Short	1 = Small
216 6072 = Thigh high plain band w/dots		04 = Java/Nightshade			3 = Long	2 = Medium
216 1002 = Pantyhose		11 = Beige				3 = Large
		17 = Black				4 = Xlarge
		22 = Espresso				

VenoTrain® discrétion THE TRANSPARENT **STOCKING**

The exceptionally fine material of these compression stockings gives them a stylish transparency and lightness.

Available in various styles and versatile colors, these stockings will give you an extra sense of vitality and well-being in your work, leisure, and everyday life.

C1

COLORS



FEATURES

- · Highly elastic, making the stockings easy to put on
- Hugs the legs without constricting them
- Breathable fine knit for a long-lasting wearing comfort

DETAILS

COMPRESSION CLASS

• Ccl 1 = 20-30 mmHg**TOE** closed **MATERIAL CONTENT**

· 70 % polyamide, 30 % elastane


CARE INSTRUCTIONS



STANDARD

	AD	AG	AT
Clair	•	•	
Beige	•	•	•
Espresso	•	•	•
Black	•	•	•
Java	•	•	•

SIZE CHART





	Ccl	Color	Toe	Circumference	Length	Size
217 800 = Knee high	1 = 20 – 30 mmHg	12 = Caramel	00 = Closed toe	2 = Normal	2 = Short	1 = Small
217 606 = Thigh high lace band	2 = 30 – 40 mmHg	15 = Navy		3 = Plus	3 = Long	2 = Medium
217 100 = Pantyhose		17 = Black				3 = Large
		69 = Anthracite				4 = Xlarge

VenoTrain® lookATTRACTIVE AND SHIMMERING

The VenoTrain look is a real highlight among medical compression stockings thanks to its glossy look – for a stylish appearance on any occasion!

The VenoTrain look enhances the figure with a push-up effect. Stretch comfort and its special trilobal fiber guarantee a perfect fit.

STANDARD COLORS



FEATURES

- Fashionable and elegant shimmer effect
- · Slightly transparent

DETAILS

COMPRESSION CLASSES

- Ccl 1 = 20-30 mmHg
- · Ccl 2 = 30-40 mmHg

TOE closed

MATERIAL CONTENT

- · Ccl 1: 60 % polyamide, 40 % elastane
- · Ccl 2: 55 % polyamide, 45 % elastane

CARE INSTRUCTIONS




STANDARD







SIZE CHART



C0 C1 C2 C3 C4 C5 C6



	Ccl	Color	Toe	Circumference	Length	Size
218 800 = Knee high	1 = 20-30 mmHg	11 = Cream	00 = Closed toe	2 = Normal	2 = Short	1 = Small
218 807 = Knee high w/silicone band 5 cm*	2 = 30-40 mmHg	12 = Caramel	20 = Open toe	3 = Plus	3 = Long	2 = Medium
218 605 = Thigh high lace band w/dots		14 = Amaretto				3 = Large
218 607 = Thigh high plain band w/dots**		15 = Navy				4 = Xlarge
218 100 = Pantyhose		17 = Black				
218 130 = Maternity		22 = Espresso				
		69 = Anthracite				
		46 = Wild Berry				

^{*} Available only in black and caramel.
** Not available in fashion colors.

Tango	Ccl	Color	Toe	Circumference	Length	Size
284 800 = Knee high	1 = 20-30 mmHg		00 = Closed toe	2 = Normal	2 = Short	1 = Small
284 605 = Thigh high lace band w/dots	2 = 30-40 mmHg	in all colors		3 = Plus	3 = Long	2 = Medium
284 100 = Pantyhose						3 = Large
						4 = Xlarge

Jive	Ccl	Color	Toe	Circumference	Length	Size
285 800 = Knee high	1 = 20-30 mmHg		00 = Closed toe	2 = Normal	2 = Short	1 = Small
285 605 = Thigh high lace band w/dots	2 = 30-40 mmHg	in all colors		3 = Plus	3 = Long	2 = Medium
285 100 = Pantyhose						3 = Large
						4 = Xlarge

VenoTrain® micro FASHIONABLE, MODERN, AND **VERSATILE**

With a particularly high microfiber content that makes it comfortable to wear, the VenoTrain micro is one of our most popular compression stockings.

Available in many designs and a wide selection of colors, it highlights the individual style of your customers.

STANDARD COLORS



ADDITIONAL COLORS



FEATURES

- More than 50 per cent microfiber content: breathable, skinfriendly, and easy-care
- Figure-shaping panty area (pantyhose) and comfort sole
- Opaque full-fashioned stocking quality

DETAILS

COMPRESSION CLASSES

- · Ccl 1 = 20-30 mmHg
- · Ccl 2 = 30-40 mmHq TOE closed, open

MATERIAL CONTENT

- · Ccl 1: 65 % polyamide, 35 % elastane
- · Ccl 2: 55 % polyamide, 45 % elastane

CARE INSTRUCTIONS



STANDARD FASHION COLORS



STANDARD










SIZE CHART

CO

C1









THE SECOND SKIN FOR TRENDSETTERS

A true fashion statement: the Design Edition and individually dyed VenoTrain micro Batik versions create real highlights.

With a subtle Tango or Jive pattern, or wildly textured in batik with navy or black, your customers' legs will be real eye-catchers.

JIVE



VenoTrain® microDESIGN EDITION AND BATIK


STANDARD COLORS



ADDITIONAL COLORS



STANDARD FASHION COLORS



STANDARD







61 62 63 64



		Ccl	Color	Toe	Circumference	Length	Size
Long	208 800 = Knee high	1 = 20-30 mmHg	12 = Caramel	00 = Closed toe	2 = Normal	2 = Short	1 = Small
Short	208 810 = Knee high	2 = 30-40 mmHg	15 = Navy		3 = Plus	3 = Long	2 = Medium
			69 = Anthracite				3 = Large
			17 = Black				4 = Xlarge

VenoTrain® business TIMELESSLY ELEGANT

While traveling or at work: the VenoTrain business meets every requirement for a premium knee-high stocking. The combination of cotton and Timbrelle® high-tech fiber offers top-class wearing comfort with thermal regulation. Highly elastic LYCRA® and the extra-wide supportive band ensure that the stocking is easy to put on and stays perfectly in place.

STANDARD COLORS



FEATURES

- · Innovative mix of materials
- Reinforced sole
- Timbrelle® high-tech fiber providing longlasting freshness and healthy foot conditions



DETAILS

COMPRESSION CLASSES

- Ccl 1 = 20-30 mmHg
- Ccl 2 = 30-40 mmHg

TOE closed

MATERIAL CONTENT

- · Ccl 1: 80 % polyamide, 14 % elastane, 6 % cotton
- · Ccl 2: 71% polyamide, 24 % elastane, 5 % cotton

CARE INSTRUCTIONS



STANDARD



SIZE CHART

	CIDCUMEED	CIRCUMFERENCE		5	1	М	I	L XL		(L
	CIRCOMFERI	ENCE	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
ength ID										
	cC	cm	28-33	33-38	33-38	38-43	38-43	43-48	43-48	48-5
ng over 41 cm (16")	Calf	inch	11 – 13	13-15	13-15	15-17	15-17	17-19	17-19	19-2
pe Knee high										
	сВ	cm	19-	- 22	22-	- 25	25-	-28	28-	-32
	Ankle	inch	7.5-	-8.5	8.5	-10	10-	-11	11-	12.5



TREATMENT IN PHLEBOLOGICAL SPECIAL CASES FOR MODERATE TO SEVERE VEIN PROBLEMS

Some patients with vein problems have **special requirements** that call for compression stockings with unique modes of action and properties. These requirements can be very diverse. That is why the VenoTrain special range offers compression stockings especially for patients with:

- GENTLE AND EFFECTIVE VenoTrain soft
- RESTRICTED MOBILITY VenoTrain soft S
- PARTICULARLY WEAK VEINS VenoTrain impuls
- OPEN WOUNDS VenoTrain ulcertec



CCI	Color	106	Circumference	Length	Size
1 = 20-30 mmHg	11 = Cream	00 = Closed toe	2 = Normal	2 = Short	1 = Small
2 = 30-40 mmHg	12 = Caramel	20 = Open toe	3 = Plus	3 = Long	2 = Medium
	14 = Amaretto				3 = Large
	15 = Navy				4 = Xlarge
	17 = Black				
	22 = Espresso				
	42 = Anthracite				
	1 = 20-30 mmHg	1 = 20-30 mmHg 11 = Cream 2 = 30-40 mmHg 12 = Caramel 14 = Amaretto 15 = Navy 17 = Black 22 = Espresso	1 = 20-30 mmHg 11 = Cream 00 = Closed toe 2 = 30-40 mmHg 12 = Caramel 20 = Open toe 14 = Amaretto 15 = Navy 17 = Black 22 = Espresso	1 = 20-30 mmHg 11 = Cream 00 = Closed toe 2 = Normal 2 = 30-40 mmHg 12 = Caramel 20 = Open toe 3 = Plus 14 = Amaretto 15 = Navy 17 = Black 22 = Espresso	1 = 20-30 mmHg 11 = Cream 00 = Closed toe 2 = Normal 2 = Short 2 = 30-40 mmHg 12 = Caramel 20 = Open toe 3 = Plus 3 = Long 14 = Amaretto 15 = Navy 17 = Black 22 = Espresso

VenoTrain® soft GENTLE AND EFFECTIVE

The VenoTrain soft, always the right choice: this versatile stocking can withstand plenty of wear, yet it is very soft and gentle on the skin.

It is available in many styles with a proven medical effect and a beneficial massage action.

STANDARD COLORS



ADDITIONAL COLORS



FEATURES

- For men and women
- · Gentle massage action
- Particularly gentle on the skin

DETAILS

COMPRESSION CLASSES

- Ccl 1 = 20-30 mmHg
- Ccl 2 = 30-40 mmHg

TOE closed, open MATERIAL CONTENT

- · Ccl 1: 70 % polyamide, 30 % elastane
- · Ccl 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS



STANDARD










SIZE CHART



CO C1 C2 C3 C4 C5



Soft arm sleeve	Ccl	Color	Length	Size
212 554 = silicone band	1 = 20-30 mmHg	00 = Natural	2 = Short	1 = Small
	2 = 30-40 mmHg	12 = Caramel	3 = Long	2 = Medium
		17 = Black		3 = Large
				4 = Xlarge

VenoTrain® soft Arm EVERYDAY MOBILITY

Seamless with an extra-soft inner elbow area. This facilitates everyday mobility with complete medical efficacy.

COLORS







FEATURES

- · Comfort elbow for greater freedom of movement, wearing comfort, and
- a reduced pressure sensation
- Noticeable effect throughout the entire wearing period
- · Gentle on the skin, soft, and breathable

DETAILS

- **MATERIAL CONTENT**
- · Ccl 1: 74 % polyamide, 26 % elastane
- · Ccl 2: 69 % polyamide, 31 % elastane

COMPRESSION CLASSES

- Ccl 1 = 20-30 mmHg
- · Ccl 2 = 30-40 mmHg


CARE INSTRUCTIONS



STANDARD



SIZE CHART



		S	М	L	XL
		min. – max.	min. – max.	min. – max.	min. – max.
cG	cm	24 - 30	27 – 33	29 – 36	32 – 40
	inch	9 - 12	11 – 13	11 – 14	13 – 16
cE	cm	22 – 26	24 – 29	27 – 32	30 – 35
	inch	9 – 10	9 – 11	11 – 13	12 – 14
cC	cm	14 – 15	15 – 16.5	16.5 – 18	18 – 20
	inch	5.5 – 6	6 – 6.5	6.5 – 7	7 – 8



	Ccl	Color	Toe	Circumference	Length	Size
212 800 = Knee high	1 = 20-30 mmHg	12 = Caramel	04 = Closed toe	2 = Normal	2 = Short	1 = Small
212 804 = Knee high w/silicone band (3 cm)	2 = 30-40 mmHg		24 = Open toe	3 = Plus	3 = Long	2 = Medium
212 607 = Thigh high plain band w/dots						3 = Large
						4 = Xlarge

VenoTrain® soft SCOMFORT IN THE EVENT OF RESTRICTED MOBILITY

Based on proven VenoTrain soft quality which offers a gentle massage effect to stimulate the blood flow, the VenoTrain soft S also takes into account special requirements of people who spend most of their time seated.

The wide heel and comfortable instep area ensure that the foot, which is angled during sitting, has enough space, and that no creases can develop.

The robust knitted fabric can withstand plenty of wear, and can also be easily put on by care staff.

STANDARD COLORS



FEATURES

- Wide heel area
- Adapted to the typical anatomical foot posture while sitting (90-degree angle)
- Breathable material with a revitalizing effect

DETAILS

COMPRESSION CLASSES

- · Ccl 1 = 20-30 mmHg
- · Ccl 2 = 30-40 mmHg

TOE closed, open MATERIAL CONTENT

- · Ccl 1: 70 % polyamide, 30 % elastane
- · Ccl 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS




STANDARD






SIZE CHART

CO





STANDARD SIZE CHART



	Ccl	Color	Toe	Circumference	Length	Size
207 800 = Knee high	2 = 30-40 mmHg	12 = Caramel	00 = Closed toe	2 = Normal	2 = Short	1
207 804 = Knee high w/silicone band (3 cm)	3 = 40-50 mmHg	17 = Black	20 = Open toe	3 = Plus	3 = Long	2
207 607 = Thigh high plain band w/dots						3
207 620 = Thigh high w/waist attachment R		-	-			4
207 630 = Thigh high w/waist attachment L						5
						6

VenoTrain® impuls POWERFUL SUPPORT FOR PARTICULARLY WEAK VEINS

Produces a particularly deep action by exerting high working pressure with even the slightest muscle movements. Easy to put on and breathable thanks to microfibers.

Highly effective for advanced vein problems, post-thrombotic syndrome, follow-up care for healed ulcers, and therapy in the early stages of lymphedema.

STANDARD COLORS



FEATURES

- The short-stretch knitted fabric generates a high working pressure with even the slightest muscle movements
- Skin-friendly microfiber material makes the product breathable and ensures temperature control

DETAILS

COMPRESSION CLASSES

- \cdot Ccl 2 = 30-40 mmHg
- Ccl 3 = 40-50 mmHg

TOE closed, open MATERIAL CONTENT

- · Ccl 2: 80 % polyamide, 20 % elastane
- · Ccl 3: 75 % polyamide, 25 % elastane

CARE INSTRUCTIONS



STANDARD










62 63 64



STANDARD SIZE CHART



	Ccl	Color	Toe	Circumference	Length	Size
209 800 = Knee high w/liners	0 = 30 – 40 mmHg	00 = Natural	21 = Open toe (w/closed toe liners)	2 = Normal	2 = Short	0 = Xsmall
				3 = Plus	3 = Long	1 = Small
						2 = Medium
209 800 = Liners	0	10	01			3 = Large
						4 = Xlarge

VenoTrain® ulcertec FOR VENOUS LEG ULCERS

The combination of an overstocking with a patented rhomboid knitted fabric and a skin-friendly liner for holding the wound dressing in place makes this compression stocking easy to put on without any need for assistance.

The liners are also ideal for providing a continuous low level of compression during the night. These are available to order separately in a pack of three and are suitable for washing at high temperatures.

COLORS



FEATURES

- · Coordinated compression of overstocking and liner
- Unique compression profile that accelerates the healing process while ensuring exceptional patient comfort
- Differentiated treatment based on two designs within one compression class
- · Liner with gentle compression providing a compression pressure in the lower range of compression class 3 in combination with the overstocking

DETAILS

COMPRESSION CLASSES

· Ccl 2 = 30-40 mmHg

TOE

- Open (overstocking)
- · Closed (liner)

MATERIAL CONTENT

· 87 % polyamide, 13 % elastane

CARE INSTRUCTIONS

· Overstocking: 40 💥









STANDARD





LYMPHOLOGICAL TREATMENT

Bauerfeind has a wide range of compression products for lymphological conditions which take into account the patients' different requirements.



Thanks to our custom-made **VenoTrain curaflow** flat knit products, we can offer lymphatic patients perfectly tailored solutions with more than 1,000 different options of combining products and additional features.

The beneficial therapeutic features of the range include the knitted fabric's high microfiber content, the flat sensitive stitching, and special comfort zones such as the extra-soft area behind the knee.



ACCESSORIES

THE PERFECT ADDITION FOR EVERYDAY USE

These handy extras from the VenoTrain accessories range are useful for ensuring that your patients feel comfortable using and wearing their VenoTrain compression stockings. They provide effective treatment support, make it easier to put on the products and care for the compression stockings.

01 Special detergent

cleans gently, ensuring careful treatment of your products, and providing lasting protection

02 Rubber gloves

help with putting the compression stockings on correctly





VenoTrain® glider COMPACT AID FOR PUTTING ON AND TAKING OFF STOCKINGS

The VenoTrain glider is a practical aid for putting on and taking off medical compression stockings. It can be prescribed by a physician. The lightweight, hard-wearing material makes the stocking glide effortlessly over the arm or leg. It is suitable for compression stockings with an open or closed toe. Saving space and time, it is the ideal companion when you're out and about.

FEATURES

- Extremely lightweight, robust, and tear-resistant
- · Easy to handle
- · Ideal for travel



PUTTING ON THE STOCKINGS: STEP 01



PUTTING ON THE STOCKINGS: STEP 0:



PUTTING ON THE STOCKINGS: STEP 03



PUTTING ON THE STOCKINGS: STEP 04



PUTTING ON THE STOCKINGS: STEP 05



TAKING OFF THE STOCKINGS

VenoTrain® glider plusROBUST DONNING AID

The VenoTrain glider plus is a robust donning aid with an ergonomic design for all types of medical compression stockings. It can be prescribed by the physician.

Folding handles make it easy to put on the stockings, avoiding the need to bend down. It has been tested by an independent institute for ease of use, even for customers with restricted mobility.

FEATURES

· A compact, ergonomic design, light metal frame with folding handles, easy to use without the need for any great effort





PUTTING ON THE STOCKINGS: STEP 04



PUTTING ON THE STOCKINGS



PUTTING ON THE STOCKINGS: STEP 05



PUTTING ON THE STOCKINGS: STEP 03



PUTTING ON THE STOCKINGS: STEP 06

A BETTER QUALITY OF LIFE FOR YOUR CUSTOMERS

Have you got any questions, suggestions, or requests? Please contact us:

GALIEN CANADA

2017 Cunard H7S 2N1

Laval, Quebec

P Office local +1 (450) 629-9889 P Office toll free +1 877 629-9889 F Local +1 (450) 629-9895 F Toll free +1 877 629-9895

E info@galiengroup.com



